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- [1] Stenman UH, Tiitinen A, Alifan H, et al. The classification, functions and clinical use of different isoforms of HCG. Hum Reprod Update, 2006, 12:769-784.
- [2] Cole LA, Sutton JM. Selecting an appropriate hCG test for managing gestational trophoblastic disease and cancer. J Reprod Med, 2004, 49:545-553.
- [3] Seki K, Matsui H, Sekiya S. Advances in the clinical laboratory detection of gestational trophoblastic disease. Clin Chim Acta, 2004, 349:1-13.
- [4] McChesney R, Wilcox AJ, O'Connor JF, et al. Intact HCG, free HCG beta subunit and HCG beta core fragment: longitudinal patterns in urine during early pregnancy. Hum Reprod, 2005, 20: 928-935.
- [5] Butler SA, Iles RK. The free monomeric beta subunit of human chorionic gonadotrophin (hCG beta) and the recently identified homodimeric beta-beta subunit (hCG beta beta) both have autocrine growth effects. Tumour Biol, 2004, 25:18-23.
- [6] Behtash N, Ghaemmaghami F, Honar H, et al. Is normal beta-hCG regression curve helpful in the diagnosis of persistent trophoblastic disease? Int J Gynecol Cancer, 2004, 14:980-983.
- [7] 王小平,向阳,张德永,等. 化学发光法测定绒毛膜促性腺激素游离  $\beta$  亚单位对滋养细胞疾病的诊断价值. 中国实用妇科与产科杂志,2001,17:660-662.

- [8] Cole LA, Khanlian SA, Muller CY, et al. Gestational trophoblastic diseases; 3. Human chorionic gonadotropin-free beta-subunit, a reliable marker of placental site trophoblastic tumors. *Gynecol Oncol*, 2006, 102:160-164.
- [9] Kohorn EI, Cole L. Nicked human chorionic gonadotropin in trophoblastic disease. *Int J Gynecol Cancer*, 2000, 10:330-335.
- [10] Zavadil M, Feyereis J, Safur P, et al. Syndrome of persistent low levels of human chorionic gonadotropin (hCG). *Ceska Gynekol*, 2006, 71:136-142.
- [11] Cole LA. Phantom hCG and phantom choriocarcinoma. *Gynecol Oncol*, 1998, 71:325-329.
- [12] Ward G, McKinnon L, Badrick T, et al. Heterophilic antibodies remain a problem for the immunoassay laboratory. *Am J Clin Pathol*, 1997, 108:417-421.
- [13] Esfandiari N, Goldberg JM. Heterophile antibody blocking agent to confirm false positive serum human chorionic gonadotropin assay. *Obstet Gynecol*, 2003, 101:1144-1146.
- [14] Olsen TG, Hubert PR, Nycum LR. Falsely elevated human chorionic gonadotropin leading to unnecessary therapy. *Obstet Gynecol*, 2001, 98:843-845.
- [15] Cole LA, Butler S. Detection of hCG in trophoblastic disease. The USA hCG reference service experience. *J Reprod Med*, 2002, 47:433-444.
- [16] Cole LA, Khanlian SA. Hyperglycosylated hCG; a variant with separate biological functions to regular hCG. *Mol Cell Endocrinol*, 2007, 260:228-236.
- [17] Cole LA, Khanlian SA, Riley JM, et al. Hyperglycosylated hCG in gestational implantation and in choriocarcinoma and testicular germ cell malignancy tumorigenesis. *J Reprod Med*, 2006, 51:919-929.
- [18] Cole LA, Butler SA, Khanlian SA, et al. Gestational trophoblastic diseases; 2. Hyperglycosylated hCG as a reliable marker of active neoplasia. *Gynecol Oncol*, 2006, 102:151-159.
- [19] Khanlian SA, Cole LA. Management of gestational trophoblastic disease and other cases with low serum levels of human chorionic gonadotropin. *J Reprod Med*, 2006, 51:812-818.

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